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This device is intended for a specialized procedure, which must be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques and always apply them under appropriate conditions, including in an operating room.

DESCRIPTION

Four exclusive Prosthetic Abutments are required for performing the One Step Hybrid Technique.

Brass Coping: a Prosthetic Abutment for laboratory use that presents a tapered geometry in the outer part for fitting the Castable Coping. In its inner part, there is a slot compatible with the planned Prosthetic Abutment.

Castable Coping: a cylindrical Prosthetic Abutment with retentions on its outer part for the waxing. In its inner part, there is a slot compatible with the planned Prosthetic Abutment.

Titanium Coping: a Prosthetic Abutment with outer retentions and treated surface (either porous or rugged) for retaining cement. It is the prosthetic intermediate of the One Step Hybrid technique that is cemented to the metallic infrastructure and screwed between the intermediate and the prosthesis.

Working Screw - Long (laboratory room): with a 1.2 mm hexagon fit in one end and thread compatible with the prosthetic intermediate diameter in the other end.

INDICATIONS FOR USE

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with appropriate occlusal loading. Multiple tooth applications may be rigidly splinted.

APPLICATIONS

The One Step Hybrid Technique is indicated for multiple screwed prostheses when there is a need to obtain passivity in the metallic infrastructure fitting without using weld and regardless of the metal dimensional change after the casting process.

They are available according to the table below:

Compatible component	Castable Coping	Brass Coping	Titanium Coping
Mini Conical Abutment Ø 4.1	118.083 Mini Conical Abutment One Step Hybrid Coping 4.1	118.081 Mini Conical Abutment One Step Hybrid Coping 4.1	118.082 Mini Conical Abutment One Step Hybrid Coping 4.1
	118.340 Neo Burn-Out Coping Mini Conical Abutment 4.1	118.331 Neo Mini Conical Abutment One Step Hybrid Coping 4.1	118.330 Neo Mini Conical Abutment One Step Hybrid Coping 4.1
Mini Conical Abutment Ø 5.0	118.089 Castable Cylinder do Mini Conical Abutment 5.0	118.087 Cylinder in Tin do Mini Conical Abutment 5.0	118.088 Cylinder in Titanium do Mini Conical Abutment 5.0
Micro Conical Abutment	118.250 Micro Conical Abutment One Step Hybrid Coping	118.248 Micro Conical Abutment One Step Hybrid Coping	118.249 Micro Conical Abutment One Step Hybrid Coping
	118.341 Neo Burn-Out Coping Micro Conical Abutment	118.333 Neo Micro Conical Abutment One Step Hybrid Coping	118.332 Neo Micro Conical Abutment One Step Hybrid Coping
Transepithelial Abutment	118.086 Transepithelial Abutment One Step Hybrid Coping	118.084 Transepithelial Abutment One Step Hybrid Coping	118.085 Transepithelial Abutment One Step Hybrid Coping
CM Abutment / GT Implant	118.174 CM Abutment One Step Hybrid	118.173 CM Abutment One Step Hybrid	118.172 CM Abutment One Step Hybrid

Check the compatibility between the chosen interfaces. The Prosthetic Abutments for the One Step Hybrid Technique are compatible with the implants/abutments from the same lines, according to previous tables.

While selecting the prosthetic abutment, pay attention to its use indication.

CONTRAINDICATIONS

Be aware of patients who present signs of allergic reaction or hyper-sensitivity to the chemical elements of the material: titanium.

This product is contraindicated for insufficient interocclusal space and unsatisfactory three-dimensional position of the implant.

HANDLING

Whenever the Prosthetic Components for the One Step Hybrid Technique are used in procedures of two stages, a previous preparation of soft tissues with a Healing Abutment is recommended.

After installing the implants/intermediates in the mouth, transfer their position through molding with their corresponding Impression Coping, according to appropriate techniques.

CAST METALLIC STRUCTURE: Place the Brass Copings over Analogs that correspond to the intermediate prosthetic abutment. On top of the Brass Copings, place the Castable Copings, fastening them with Working Screws (long). After that, connect with acrylic resin the Castable Copings that are parallelly connected.

Perform the infrastructure waxing and, once it is concluded, remove the Working Screws, submitting the infrastructure to the casting process with specific metallic alloy. Place the casting infrastructure over the Brass Copings in order to check the fitting passivity. If necessary, provide internal wear in the regions corresponding to the Castable Copings, which are now casted in metal, aiming to achieve infrastructure passivity on Brass Copings. Once the infrastructure passivity is achieved, perform internal retentions in the casted Copings in order to create mechanical retentions for the cement. Apply a thin layer of primer specific for metal (alloy primer) on this area. Substitute the Brass Copings, which present a slightly greater dimension (0.10 mm), for Titanium Copings, fastening them over their corresponding Analogs with their corresponding Screws (short). Apply the primer specific for metal (alloy primer) on the external portion of Titanium Copings, obliterating the input hole with wax to avoid the entry of resinous cement. Also on the outer surface of the Titanium Copings and on the internal portion that corresponds to the Coping in the infrastructure, apply a dual resinous cement. Still with fresh cement, press the infrastructure over the Titanium Copings, immediately removing any overflow excess from the orifice. After cementing it, unscrew the infrastructure from the model and remove all the excess of remaining cement from the Titanium Copings edge. **NOTE:** For the cementation, it is recommended to use dual resinous cement Panavia F (Kuraray Co Ltd Tokyo-Japan) and the primer specific for metal Alloy Primer (Kuraray Co Ltd Tokyo-Japan). Proof and tests of passivity and adaptation of the prosthesis structure must be performed. The connections used for handling of Prosthetic Abutments for the One Step Hybrid Technique must follow the table below:

Prosthetic Abutments for the One Step Hybrid Technique	Recommended Connection	
118.081 Mini Conical Abutment One Step Hybrid Coping 4.1	1.2 Connection	
118.082 Mini Conical Abutment One Step Hybrid Coping 4.1		
118.082-ET Mini Conical Abutment One Step Hybrid Coping 4.1		
118.083 Mini Conical Abutment One Step Hybrid Coping 4.1		
118.084 Transepithelial Abutment One Step Hybrid Coping		
118.085 Transepithelial Abutment One Step Hybrid Coping		
118.086 Transepithelial Abutment One Step Hybrid Coping		
118.087 Mini Conical Abutment One Step Hybrid Coping 5.0		
118.088 Mini Conical Abutment One Step Hybrid Coping 5.0		
118.088-ET Mini Conical Abutment One Step Hybrid Coping 5.0		
118.089 Mini Conical Abutment One Step Hybrid Coping 5.0		
118.172 CM Abutment One Step Hybrid		
118.173 CM Abutment One Step Hybrid		
118.174 CM Abutment One Step Hybrid		
118.248 Micro Conical Abutment One Step Hybrid Coping		
118.249 Micro Conical Abutment One Step Hybrid Coping		
118.249-ET Micro Conical Abutment One Step Hybrid Coping		
118.250 Micro Conical Abutment One Step Hybrid Coping		
118.330 Neo Mini Conical Abutment One Step Hybrid Coping 4.1		Neo Connection
118.330-ET Neo Mini Conical Abutment One Step Hybrid Coping 4.1		
118.331 Neo Mini Conical Abutment One Step Hybrid Coping 4.1		
118.332 Neo Micro Conical Abutment One Step Hybrid Coping		
118.332-ET Neo Micro Conical Abutment One Step Hybrid Coping		
118.333 Neo Micro Conical Abutment One Step Hybrid Coping		
118.340 Neo Burn-Out Coping Mini Conical Abutment 4.1		
118.341 Neo Burn-Out Coping Micro Conical Abutment		

The torque recommended for the screws of all Copings is of 10 N.cm.

PRESENTATION AND STERILIZATION

These products are for single use only and supplied non-sterile, being unitarily packaged.

PRECAUTIONS

- It is recommended that the final prosthesis is cleaned before installation in the mouth, according to laboratory instructions.
- Before installing the product, make sure that it has the indication for use of the selected prosthetic abutment. Make sure that the stability of the implant is sufficient to support the prosthetic abutment installation torque and functional load, in accordance with the instructions for use of the implant.
- The material selection of the prosthesis structure must consider general aspects of the patient.
- Do not use the product after its expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic abutment used. The excess or lack of torque can lead to undesirable results.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- During installation, make sure to align it with the insertion axis of the implant. Make sure it is perfectly seated on the implant. To do so, it is recommended to use periapical X-rays with the parallel technique.
- Make sure to use a prosthetic screw that is compatible with the prosthetic interface and with the prosthetic abutment.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, deleted markings, compromised sharpening, deformation or wear.
- It is the dental surgeon's responsibility to use the Neodent products according to the instructions for use.
- Do not use the product if the packaging is violated.
- This product is for single use only and cannot be reused.
- Reuse of this product may cause: • adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing; • changes in physical,

mechanical and chemical properties of products, macro and micro structural, that can put at risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.

- The result shall only be achieved if the sequence of Neodent products is used. The utilization of instruments and/or implants and/or prosthetic abutments from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product warranty.
- The improper planning may endanger the implant/prosthesis set performance resulting in system faults, such as Implant loss or fracturing, loosening or fracturing of the Prosthetic Screws and/or metallic infrastructure, or releasing of the titanium Coping from the metallic infrastructure.

ADVERSE EFFECTS

Adverse effects are not expected as long as the product is used according to its instructions for use.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need for a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. The professional in charge is responsible for providing these guidelines.

STORAGE CONDITIONS

This product must be stored, in its original packaging, in a clean and dry location, at a maximum temperature of 40 °C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them. Before discarding them into the environment, it is recommended to observe the current legislation and adhere to it.


DATE OF EXPIRATION

Written on the label.

SYMBOLS

The table describes the symbols that may be printed on the product labeling. Please, refer to the physical labeling to see which symbols are applicable to the product.

Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture

Symbols	Description
	Manufacturer
	Use-by date
Qty	Quantity
	Consult instructions for use
PROIBIDO REPROCESSAR	Reprocessing forbidden
STERILE EO	Sterilized using ethylene oxide
STERILE R	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
EC REP	Authorized representative in the European Community
EC IMP	Importer in the European Community
CH REP	Swiss Authorised representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
Rx only	Notification required by FDA for United States market
	CE Mark
NON-STERILE 	Non-sterile

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Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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© 2021 - JJGC Indústria e Comércio de Materiais Dentários S.A. All rights reserved.
© 2021 - JJGC Indústria e Comércio de Materiais Dentários S.A. Med ensamrätt.