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**INDICATIONS FOR USE**

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore single-stage or two-stage surgical procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with physiological occlusal loading. Multiple teeth applications can be rigidly splinted.

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According to the image below, the transfer piece presents several features: (1) an hexagon compatible with the Neodent Hexagonal Implant Driver – this feature is used to apply torque during the installation; (2) an intermediate safety feature (fuse) – intended to break if a torque higher than 60 N.cm is applied to the implant in order to protect the integrity of the implant; (3) a secondary hexagon for implant removal (counter-clockwise).

**APPLICATIONS**

The Zirconia Implant is indicated for surgical intraoral installation in bone with density I / II or in bone III / IV, according to Lekholm & Zarb's bone classification (1985), in the maxilla or mandible. It can be used as a support for single-unit or multiple-unit prostheses in immediate or conventional loading protocol. It may be installed immediately after the extraction of dental root. **Note:** For the purposes of immediate loading, primary stability must reach, at least, 35 N.cm and the patient must present physiological occlusion.

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**WARNING**

The maximum insertion torque for the Zirconia Implant is 60 N.cm. Applying a torque higher than 60 N.cm may cause damages and/or break the implant. If maximum torque is reached and the insertion cannot be concluded, it is recommended to remove the implant and reprepare the implant bed for a new insertion attempt.

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### PRECAUTIONS

- Due to limited mouth opening in the posterior region, an assessment is required prior to installation of long implants in the region of premolars and molars. The installation of long implants in these regions may require them to be installed at an angle.
- Inadequate surgical and/or prosthetic planning can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of abutments and/or prosthetic screws.
- Be aware in cases of patients that present signs of allergy or hypersensitivity to chemical elements of the material: stainless steel.
- Observe the conditions of the intraoral tissue, bone quality and quantity of the bed receiving the implant, through radiographic and/or tomography examinations. Not performing the pre-surgical assessment may compromise the success of the procedure.
- As for the systemic aspect, consider the general health of the patient. In particular, one must be careful in cases of patients who have allergies to drugs, local or systemic factors that may interfere with the healing process of the bone tissues or soft tissues, or the process of osseointegration. For example, the bone already exposed to radiation in the head and neck area, diabetes mellitus, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, previous treatment of jaws pathologies and oral mucosa abnormalities.
- Treatment with bisphosphonates results in potential risk of peri-implant osteonecrosis.
- Avoid approaching the proximities of the mandibular nerve channel during bone bed surgical preparation. Nerve damage may result in anesthesia, paresthesia, and dyesthesia.
- The material to be used during the procedure must be sterile. Never use potentially contaminated items, since contamination may lead to infections.

- This product is for single use and cannot be re-sterilized. Neodent does not accept any responsibility for re-sterilized implants, regardless of who has carried out or how it was performed. Cleaning, disinfection, and sterilization of this product by oneself may compromise essential material and design characteristics.
- Reuse of this product may cause: adverse biological effects of residual products, microorganisms, and/or substances resulting from previous uses and/or reprocessing; changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put at risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if its packaging is damaged.
- Do not use the product with its validity expired.
- The sterilization is guaranteed only if the sterile barrier (full blister) is not damaged.
- This product must be used immediately after opening of its packaging, at the moment of surgery. If it is not used, discard it.
- Follow the rules of asepsis when removing the implant from the sterile packaging.
- The maximum installation torque suggested is of 60 N.cm. Insertion torque greater than the one recommended can make the system inoperative. Applying an improper torque could lead to implant breakage or bone necrosis.
- Before each procedure, make sure the pieces are properly seated. Ensure the perfect fit between implant and transfer piece.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Do not continue to insert the implant in case the transfer piece breaks. Insisting on the installation may cause damages to the system and/or the patient. Use the appropriate connection to remove the implant.
- Before each procedure, check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is any damage, deleted markings, compromised sharpening, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the dental surgeon's responsibility to use the Neodent products according to the instructions for use.
- In case there is an accident related to the device, inform the manufacturer using the communication channels available on the website.

**ADVERSE EFFECTS**

The installation of dental implants, as well as any other surgical procedure, may cause a slight discomfort and localized edema. More persistent symptoms can occur such as: chronic pain related to the dental implant, permanent paresthesia, dysesthesia, loss of maxillary/mandibular bone edge, systemic or localized infection, oroantral or oronasal fistula, adjacent teeth unfavorably affected, irreversible damage to adjacent teeth,

fracturing of the implant, jaw, bone or prosthesis, esthetic problems, injury of the nerves, exfoliation, hyperplasia.

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- Static magnetic field of 1.5T and 3T, only.
  - Maximum spatial gradient magnetic field of 4,000 gauss/cm (40 T/m).
  - Maximum MR system indication: whole body averaged specific absorption rate (SAR) of 2 W/kg and head average SAR of 3.2 W/kg for 15 minutes of scanning (i.e. per pulse sequence) in the Normal Operating Mode.
- Under the scan conditions defined, the implants, as well as clinically relevant implant constructs, from the Neodent Implant System are expected to produce a maximum temperature rise of 4.9°C after 15-minutes of continuous scanning (i.e., per pulse sequence).
- In non-clinical testing, the image artifact caused by the implants from the Neodent Implant System extend approximately 10 mm from this device when imaged with a gradient echo pulse sequence and a 3T MR system.

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Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
	Quantity
	Consult instructions for use
	Reprocessing forbidden
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community
	Importer in the European Community
	Swiss Authorised representative
	Single sterile barrier system

Symbols	Description
	Single sterile barrier system with protective packaging inside
Rx only	Notification required by FDA for United States market
	CE Mark
NON-STERILE 	Non-sterile

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**WARNING**

Small diameter implants with angled abutments are not recommended for the molar region. The maximum insertion torque for the Zirconia Implant is 60 N.cm. Applying a torque higher than 60 N.cm may cause damages and/or break the implant. If maximum torque is reached and the insertion cannot be concluded, it is recommended to remove the implant and reprepare the implant bed for a new insertion attempt. If one insists to install the implant with a higher torque, the transfer piece has a safety feature that will fracture (2). This fracture impedes the user to continue with the implant insertion. From this moment on, implant integrity is no longer guaranteed and Neodent does not take responsibility for damages caused by the incorrect usage of the item. The removal of the implant is mandatory, which must be performed through the secondary hexagonal fitting in a counter-clockwise movement (3). After the removal, the substitution for a new implant is indicated. The secondary hexagon shall only be used for implant removal. Never use this feature to insert the implant deeper. Persisting on the insertion may cause damages to the implant and the user.

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**HANDLING**

Perform drilling of the bone bed using drills in good cutting conditions and rotation between 500 and 800 rpm for bone type III / IV, and rotation between 800 and 1200 rpm for bone type I / II, under abundant irrigation. Select the sequence of drills according to the selected implant. The insertion depth of the drills must comply with the planning of the implant final position. For the installation of Zirconia Implants, perform it manually with the Torque Wrench Driver or perform it with surgical motor using the maximum torque of 35 N.cm and drilling rotation of 30 rpm. Neodent Drills/ Bone Taps must be used in accordance with their own instructions for use. For more information, consult the Neodent catalog.  
**Sequence for Zirconia Implants:**

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**NOTES:**


The Initial Drill and Drills 2.0, 3.5, 3.75 and 4.3 indicated in the drilling sequence can be from GM or Alvim lines. Drill the bone bed according to the drill sequence indicated for the planned implant. The Bone Tap is recommended to finish the preparation of the surgical cavity in type I / II bone and should be used according to the following guidelines:

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- The cardboard and blister packagings must be opened, manually, without the use of sterile gloves.
- Break the seal of the cardboard packaging and remove the blister.
- Open the blister pack.
- Deposit the sterile flask over the surgical field. **NOTE:** The clear tube and implant must be handled with a sterile surgical glove, in a surgical environment.
- Hold the bottle using the non-dominant hand and take the lid off.
- The internal support containing the implant and transfer piece must come out attached to the lid. To do so, remove the lid and the clear tube's internal support in the axial direction without making any lateral movements.
- Using the non-dominant hand, press the sides of the internal support promoting a "pincer effect" and immobilizing the implant.
- Keep the support pressed and remove the lid.
- For installation, capture the implant transfer piece with the Hexagonal Connection, keeping it stable and slightly rotating the internal support, searching for the perfect fit between connection and transfer piece.
- Take the transfer-implant assembly to the surgical cavity.
- In the surgical motor, use maximum torque of 35 N.cm and rotation of 30 rpm.
- Complete the installation of the implant to the bone level with the Torque Wrench Driver with the Hexagonal Connection for Torque Wrench, short or long.
- At the end of the installation, make sure that one of the six hexagonal faces of the transfer piece, corresponding to the implant indexer, is facing the vestibular face.
- After the installation, remove the transfer piece from the implant.
- Load application according to torque is described in the following table:

Load Application	Min. Torque (N.cm)	Max. Torque (N.cm)
Delayed Loading*	10	60
Immediate Loading	35	60

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The model, diameter, length, position, and quantity of implants must be selected for each clinical case, considering the anatomy, the region to be rehabilitated, bone quality and quantity, as well as available space. Execute the wax-up/diagnostic of the clinical case, ensuring the appropriate alignment between implant(s), prosthetic abutment(s), and prosthesis(es).

### PRECAUTIONS

- The installation of implants with a diameter of 3.75 mm and a length of 10 mm is recommended only in the region of incisors.
- Due to limited mouth opening in the posterior region, an assessment is required prior to installation of long implants in the region of premolars and molars. The installation of long implants in these regions may require them to be installed at an angle.
- Inadequate surgical and/or prosthetic planning can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of abutments and/or prosthetic screws.
- Be aware in cases of patients that present signs of allergy or hypersensitivity to chemical elements of the material: stainless steel.
- Observe the conditions of the intraoral tissue, bone quality and quantity of the bed receiving the implant, through radiographic and/or tomography examinations. Not performing the pre-surgical assessment may compromise the success of the procedure.
- As for the systemic aspect, consider the general health of the patient. In particular, one must be careful in cases of patients who have allergies to drugs, local or systemic factors that may interfere with the healing process of the bone tissues or soft tissues, or the process of osseointegration. For example, the bone already exposed to radiation in the head and neck area, diabetes mellitus, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, previous treatment of jaws pathologies and oral mucosa abnormalities.
- Treatment with bisphosphonates results in potential risk of peri-implant osteonecrosis.
- Avoid approaching the proximities of the mandibular nerve channel during bone bed surgical preparation. Nerve damage may result in anesthesia, paresthesia, and dysthesia.

- The material to be used during the procedure must be sterile. Never use potentially contaminated items, since contamination may lead to infections.
- This product is for single use and cannot be re-sterilized. Neodent does not accept any responsibility for re-sterilized implants, regardless of who has carried out or how it was performed. Cleaning, disinfection, and sterilization of this product by oneself may compromise essential material and design characteristics.
- Reuse of this product may cause: adverse biological effects of residual products, microorganisms, and/or substances resulting from previous uses and/or reprocessing; changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put at risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if its packaging is damaged.
- Do not use the product with its validity expired.
- The sterilization is guaranteed only if the sterile barrier (full blister) is not damaged.
- This product must be used immediately after opening of its packaging, at the moment of surgery. If it is not used, discard it.
- Follow the rules of asepsis when removing the implant from the sterile packaging.
- The maximum installation torque suggested is of 60 N.cm. Insertion torque greater than the one recommended can make the system inoperative. Applying an improper torque could lead to implant breakage or bone necrosis.
- Before each procedure, make sure the pieces are properly seated. Ensure the perfect fit between implant and transfer piece.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Do not continue to insert the implant in case the transfer piece breaks. Insisting on the installation may cause damages to the system and/or the patient. Use the appropriate connection to remove the implant.
- Before each procedure, check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is any damage, deleted markings, compromised sharpening, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the dental surgeon's responsibility to use the Neodent products according to the instructions for use.
- In case there is an accident related to the device, carefully remove the parts from the patient's mouth and inform your local distributor.

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Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
	Quantity
	Consult instructions for use
	Reprocessing forbidden
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize
	Fragile

Symbols	Description
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community
	Importer in the European Community
	Swiss Authorised representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
Rx only	Notification required by FDA for United States market
	CE Mark
NON-STERILE 	Non-sterile

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This device is intended for a specialized procedure, which must be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques and always apply them under appropriate conditions, including in an operating room.

**INDICATIONS FOR USE**

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage surgical procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with physiological occlusal loading. Multiple teeth applications can be rigidly splinted.

**DESCRIPTION**

The Zirconia Implant System is composed of endosseous implants with corresponding abutments, healing abutments, cover screws, prosthetic or surgical instruments. The Zirconia Implant is manufactured in Zirconia Y-TZP (Yttria-stabilized Tetragonal Zirconia Polycrystal).  
**MACROGEOMETRY:** • Prosthetic interface with internal indexer (Zilock); • One size prosthetic platform; • Threads of trapezoidal profile; • Tapered format.  
**SURFACE:** rough surface obtained by double treatment: abrasive blasting and acid subtraction.

Implant Diameter (mm)	Implant Length (mm)
3.75	10, 11.5 and 13
4.3	

Neodent Zirconia Implants are provided with a transfer piece manufactured in stainless steel, which has the function to transfer the torque applied by the connection to the implant. The transfer piece is compatible with the Neodent Hexagonal Connection. According to the image below, the transfer piece presents several features: (1) an hexagon compatible with the Neodent Hexagonal Implant Driver – this feature is used to apply torque during the installation; (2) an intermediate safety feature (fuse) – intended to break if a torque higher than 60 N.cm is applied to the implant in order to protect the integrity of the implant; (3) a secondary hexagon for implant removal (counter-clockwise).

**APPLICATIONS**

The Zirconia Implant is indicated for surgical intraoral installation in bone with density I / II or in bone III / IV, according to Lekholm & Zarb's bone classification (1985), in the maxilla or mandible. It can be used as a support for single-unit or multiple-unit prostheses in immediate or conventional loading protocol. It may be installed immediately after the extraction of dental root. **Note:** For the purposes of immediate loading, primary stability must reach, at least, 35 N.cm and the patient must present physiological occlusion. Implants with diameter 3.75 are indicated only for the region of incisors and canines. Zirconia Implants are compatible only with prosthetic abutments from the Zirconia Implant System.

**WARNING**

Small diameter implants and angled abutments are not recommended for the posterior region. The maximum insertion torque for the Zirconia Implant is 60 N.cm. Applying a torque higher than 60 N.cm may cause damages and/or break the implant. If maximum torque is reached and the insertion cannot be concluded, it is recommended to remove the implant and reprepare the implant bed for a new insertion attempt. If one insists to install the implant with a higher torque, the transfer piece has a safety feature that will fracture (2). This fracture impedes the user to continue with the implant insertion. From this moment on, implant integrity is no longer guaranteed and Neodent does not take responsibility for damages caused by the incorrect usage of the item. The removal of the implant is mandatory, which must be performed through the secondary hexagonal fitting in a counter-clockwise movement (3). After the removal, the substitution for a new implant is indicated. The secondary hexagon shall only be used for implant removal. Never use this feature to insert the implant deeper. Persisting on the insertion may cause damages to the implant and the user.

**CONTRAINDICATIONS**

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical elements of the material: zirconia (Y-TZP), zirconium dioxide (ZrO2), yttrium oxide (Y2O3) hafnium dioxide (HfO2), and aluminum oxide (Al2O3). It is contraindicated to install the implant in the presence of acute inflammatory or infectious processes in live tissue, unsuitable bone volume and/or quality, systemic complications or diseases such as; bone metabolism disorders, blood clotting disorders, unsuitable healing capacity, incomplete jawbone growth, uncooperative and not motivated patient, abuse of drugs or alcohol, psychosis, prolonged functional disorders which resist any treatment with medications, xerostomia, weakened immunological system, diseases which require the use of steroids, endocrine diseases, insufficient oral hygiene, and pregnancy.

**HANDLING**

Perform drilling of the bone bed using drills in good cutting conditions and rotation between 500 and 800 rpm for bone type III / IV, and rotation between 800 and 1200 rpm for bone type I / II, under abundant irrigation. Select the sequence of drills according to the selected implant. The insertion depth of the drills must comply with the planning of the implant final position. For the installation of Zirconia Implants, perform it manually with the Torque Wrench Driver or perform it with surgical motor using the maximum torque of 35 N.cm and drilling rotation of 30 rpm. Neodent Drills/Bone Taps must be used in accordance with their own instructions for use. For more information, consult the Neodent catalog.

Sequence for Zirconia Implants:

Ø Implant (mm)	Bone Type	Initial Drill	Drill 2.0	Drill 3.5	Drill 3.75	Countersink Drill for Zirconia 3.75	Bone Tap for Zirconia 3.75	Drill 4.3	Countersink Drill for Zirconia 4.3	Bone Tap for Zirconia 4.3
3.75	I / II	*	✓	✓	✓	✓	✓			
	III	*	✓	✓	✓	✓				
	IV	*	✓	✓	✓					
4.3	I / II	*	✓	✓				✓	✓	✓
	III	*	✓	✓				✓	✓	
	IV	*	✓	✓				✓		

Obs: In order to prepare the surgical alveolus after extraction, use sequences of the drill used in type I bone. For the mandible, use a bone tap.  
 \* Optional  
 ✓ Recommended sequence

**NOTES:**


The Initial Drill and Drills 2.0, 3.5, 3.75 and 4.3 indicated in the drilling sequence can be from GM or Alvim lines.

Drill the bone bed according to the drill sequence indicated for the planned implant. The Bone Tap is recommended to finish the preparation of the surgical cavity in type I / II bone and should be used according to the following guidelines:

**Step 1:** To start the insertion with the Bone Tap, use the Contra-Angle handpiece. Fit the Bone Tap to the Hex Connection for Contra-Angle and configure the surgical motor with a rotation speed of a maximum of 30 rpm and a maximum torque of 35 N.cm. Start the engine and insert the Bone Tap into the surgical cavity, maintaining the drilling shaft until stability is obtained and / or it reaches the maximum torque of 35 N.cm.

**Step 2:** Then, proceed with the installation with the Torque Wrench, using the Hexagonal Connection for Torque Wrench. To do this, fit the Bone Tap to the Hexagonal Connection for Torque Wrench and perform insertion movements, clockwise, lightly pressing the Torque Wrench. If you encounter a lot of resistance, it is recommended to apply contratorque to continue with the insertion. Proceed with the insertion until the Bone Tap for Zirconia Implant reaches the mark corresponding to the chosen implant. To completely remove the Bone Tap from the surgical cavity, reverse the direction of the Torque Wrench counterclockwise and carefully remove it. If performed differently, removing the instrument can impair the thread formation.

### INSTALLATION OF THE DENTAL IMPLANT

**WARNING:** The images are merely illustrative and do not represent the actual dimensions and features of the product.

- The cardboard and blister packagings must be opened, manually, without the use of sterile gloves.
- Break the seal of the cardboard packaging and remove the blister.
- Open the blister pack.
- Deposit the sterile flask over the surgical field. NOTE: The clear tube and implant must be handled with a sterile surgical glove, in a surgical environment.
- Hold the bottle using the non-dominant hand and take the lid off.
- The internal support containing the implant and transfer piece must come out attached to the lid. To do so, remove the lid and the clear tube's internal support in the axial direction without making any lateral movements.
- Using the non-dominant hand, press the sides of the internal support promoting a "pincer effect" and immobilizing the implant.
- Keep the support pressed and remove the lid.
- For installation, capture the implant transfer piece with the Hexagonal Connection, keeping it stable and slightly rotating the internal support, searching for the perfect fit between connection and transfer piece.
- Take the transfer-implant assembly to the surgical cavity.
- In the surgical motor, use maximum torque of 35 N.cm and rotation of 30 rpm.
- Complete the installation of the implant to the bone level with the Torque Wrench Driver with the Hexagonal Connection for Torque Wrench, short or long.
- At the end of the installation, make sure that one of the six hexagonal faces of the transfer piece, corresponding to the implant indexer, is facing the vestibular face.
- After the installation, remove the transfer piece from the implant.
- Load application according to torque is described in the following table:

Load Application	Min. Torque (N.cm)	Max. Torque (N.cm)
Delayed Loading*	10	60
Immediate Loading	35	60

\* Associated with the use of Healing Abutment. When the installation torque is less than 10 N.cm, the use of the Cover Screw is recommended.

### PROSTHETIC SEQUENCE

In the immediate postoperative period, according to the chosen loading protocol, install the cover screw, healing abutment or prosthetic abutment, respecting their corresponding indications and limitations.

### TRACEABILITY LABEL

This product is supplied along with three labels that allow its traceability and must be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

### PRESENTATION AND STERILIZATION

This product is for single use and supplied sterile by the ethylene oxide method, being unitarily packed in a packaging that offers triple protection: carton box, blister and clear tube.

### PRE-OPERATORY AND PROSTHETIC PLANNING

The model, diameter, length, position, and quantity of implants must be selected for each clinical case, considering the anatomy, the region to be rehabilitated, bone quality and quantity, as well as available space. Execute the wax-up diagnostic of the clinical case, ensuring the appropriate alignment between implant(s), prosthetic abutment(s), and prosthesis(es).

### PRECAUTIONS

- The Ø3.75mm implants with angled abutments are recommended for the incisors region only.
- Due to limited mouth opening in the posterior region, an assessment is required prior to installation of long implants in the region of premolars and molars. The installation of long implants in these regions may require them to be installed at an angle.
- Inadequate surgical and/or prosthetic planning can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of abutments and/or prosthetic screws.
- Be aware in cases of patients that present signs of allergy or hypersensitivity to chemical elements of the material: stainless steel.
- Observe the conditions of the intraoral tissue, bone quality and quantity of the bed receiving the implant, through radiographic and/or tomography examinations. Not performing the pre-surgical assessment may compromise the success of the procedure.
- As for the systemic aspect, consider the general health of the patient. In particular, one must be careful in cases of patients who have allergies to drugs, local or systemic factors that may interfere with the healing process of the bone tissues or soft tissues, or the process of osseointegration. For example, the bone already exposed to radiation in the head and neck area, diabetes mellitus, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, previous treatment of jaws pathologies and oral mucosa abnormalities.
- Treatment with bisphosphonates results in potential risk of peri-implant osteonecrosis.

- Avoid approaching the proximities of the mandibular nerve channel during bone bed surgical preparation. Nerve damage may result in anesthesia, paresthesia, and dyesthesia.
- The material to be used during the procedure must be sterile. Never use potentially contaminated items, since contamination may lead to infections.
- This product is for single use and cannot be re-sterilized. Neodent does not accept any responsibility for re-sterilized implants, regardless of who has carried out or how it was performed. Cleaning, disinfection, and sterilization of this product by oneself may compromise essential material and design characteristics.
- Reuse of this product may cause: adverse biological effects of residual products, microorganisms, and/or substances resulting from previous uses and/or reprocessing; changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put at risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if its packaging is damaged.
- Do not use the product with its validity expired.
- The sterilization is guaranteed only if the sterile barrier (full blister) is not damaged.
- This product must be used immediately after opening of its packaging, at the moment of surgery. If it is not used, discard it.
- Follow the rules of asepsis when removing the implant from the sterile packaging.
- The maximum installation torque suggested is of 60 N.cm. Insertion torque greater than the one recommended can make the system inoperative. Applying an improper torque could lead to implant breakage or bone necrosis.
- Before each procedure, make sure the pieces are properly seated. Ensure the perfect fit between implant and transfer piece.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Do not continue to insert the implant in case the transfer piece breaks. Insisting on the installation may cause damages to the system and/or the patient. Use the appropriate connection to remove the implant.
- Before each procedure, check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is any damage, deleted markings, compromised sharpening, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the dental surgeon's responsibility to use the Neodent products according to the instructions for use.
- In case there is an accident related to the device, inform the manufacturer using the communication channels available on the website.

**ADVERSE EFFECTS**

The installation of dental implants, as well as any other surgical procedure, may cause a slight discomfort and localized edema. More persistent symptoms can occur such as: chronic pain related to the dental implant, permanent paresthesia, dysesthesia, loss of maxillary/mandibular bone edge, systemic or localized infection, oroantral or oronasal fistula, adjacent teeth unfavorably affected, irreversible damage to adjacent teeth, fracturing of the implant, jaw, bone or prosthesis, esthetic problems, injury of the nerves, exfoliation, hyperplasia.

Any failure in osseointegration and loss of the prosthesis during treatment may be caused by: inadequate osteotomy, infections, diseases or systemic problems, low quality or volume of the remaining bone, absence or failure of irrigation, use of nonspecific instruments and/or instruments without cutting power, insufficient oral hygiene, occlusal trauma, lack of prosthetic passivity and lack of specific training.

**MAGNETIC RESONANCE IMAGING (MRI) – SAFETY INFORMATION**

The Zirconia Implants are MRI-safe, since this is a non-magnetic ceramic material non-magnetic and, therefore, does not promote an increase in radio frequency energy and/or magnetic fields during the examination of the patient.

When the de implant is associated to abutments with a metallic screw, non-clinical testing and magnetic resonance (MRI) simulations have shown that Neodent implantable products made of metallic material are conditional to MRI (MR Conditional). A patient that has an implant from this family can be scanned safely under the following conditions:

- Static magnetic field of 1.5T and 3T, only.
- Maximum spatial gradient magnetic field of 4,000 gauss/cm (40 T/m).
- Maximum MR system indication: whole body averaged specific absorption rate (SAR) of 2 W/kg and head average SAR of 3.2 W/kg for 15 minutes of scanning (i.e. per pulse sequence) in the Normal Operating Mode.

Under the scan conditions defined, the implants, as well as clinically relevant implant constructs, from the Neodent Implant System are expected to produce a maximum temperature rise of 4.9°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implants from the Neodent Implant System extend approximately 10 mm from this device when imaged with a gradient echo pulse sequence and a 3T MR system.

**POST-OPERATIVE PRECAUTIONS AND MAINTENANCE**

Instruct the patient as to the need for professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. Evaluate the need for follow-up with periodic clinical appointments and radiographic examinations. The professional in charge is responsible for providing these guidelines.

**STORAGE CONDITIONS**

This product must be stored, in its original packaging, in a clean and dry location and protected from direct sunlight.

**DISPOSAL OF MATERIAL**

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them. Before discarding them into the environment, it is recommended to observe the current legislation and adhere to it.

**DATE OF EXPIRATION**

Written on the label.

**SYMBOLS**

The table describes the symbols that may be printed on the product labeling. Please, refer to the physical labeling to see which symbols are applicable to the product.

Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
	Quantity
	Consult instructions for use
	Reprocessing forbidden
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community
	Importer in the European Community

Symbols	Description
	Swiss Authorised representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
Rx only	Notification required by FDA for United States market
	CE Mark
NON-STERILE 	Non-sterile

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Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician. Not all products are available in all countries. Please, contact the authorized distributor.

This device is intended for a specialized procedure, which must be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques and always apply them under appropriate conditions, including in an operating room.

**INDICATIONS FOR USE**

The Neodent Implant System is intended to be surgically placed in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, to restore chewing function. It may be used with single-stage or two-stage surgical procedures, for single or multiple unit restorations, and may be loaded immediately when good primary stability is achieved and with physiological occlusal loading. Multiple teeth applications can be rigidly splinted.

**DESCRIPTION**

The Zirconia Implant System is composed of endosseous implants with corresponding abutments, healing abutments, cover screws, prosthetic or surgical instruments. The Zirconia Implant is manufactured in Zirconia Y-TZP (Yttria-stabilized Tetragonal Zirconia Polycrystal).  
**MACROGEOMETRY:** • Prosthetic interface with internal indexer (Zilock); • One size prosthetic platform; • Threads of trapezoidal profile; • Tapered format; • Bone compression capacity during installation.  
**SURFACE:** rough surface obtained by double treatment: abrasive blasting and acid subtraction.

Implant Diameter (mm)	Implant Length (mm)
3.75	10, 11.5 and 13
4.3	

Neodent Zirconia Implants are provided with a transfer piece manufactured in stainless steel, which has the function to transfer the torque applied by the connection to the implant. The transfer piece is compatible with the Neodent Hexagonal Connection.

According to the image below, the transfer piece presents several features: (1) an hexagon compatible with the Neodent Hexagonal Implant Driver – this feature is used to apply torque during the installation; (2) an intermediate safety feature (fuse) – intended to break if a torque higher than 70 N.cm is applied to the implant in order to protect the integrity of the implant; (3) a secondary hexagon for implant removal (counter-clockwise).

**APPLICATIONS**

The Zirconia Implant is indicated for surgical intraoral installation in bone with density I / II or in bone III / IV, according to Lekholm & Zarb's bone classification (1985), in the maxilla or mandible. It can be used as a support for single-unit or multiple-unit prostheses in immediate or conventional loading protocol. It may be installed immediately after the extraction of dental root. Note: For the purposes of immediate loading, primary stability must reach, at least, 35 N.cm and the patient must present physiological occlusion.

Implants with diameter 3.75 are indicated only for the region of incisors and canines. Zirconia Implants are compatible only with prosthetic abutments from the Zirconia Implant System.

**WARNING**

Small diameter implants with angled abutments are not recommended for the molar region. The maximum insertion torque for the Zirconia Implant is 60 N.cm. Applying a torque higher than 60 N.cm may cause damages and/or break the implant. If maximum torque is reached and the insertion cannot be concluded, it is recommended to remove the implant and reprepare the implant bed for a new insertion attempt.

If one insists to install the implant with a higher torque, the transfer piece has a safety feature that will fracture (2). This fracture impedes the user to continue with the implant insertion. From this moment on, implant integrity is no longer guaranteed and Neodent does not take responsibility for damages caused by the incorrect usage of the item.

The removal of the implant is mandatory, which must be performed through the secondary hexagonal fitting in a counter-clockwise movement (3). After the removal, the substitution for a new implant is indicated.

The secondary hexagon shall only be used for implant removal. Never use this feature to insert the implant deeper. Persisting on the insertion may cause damages to the implant and the user.

**CONTRAINDICATIONS**

This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical elements of the material: zirconia (Y-TZP), zirconium dioxide (ZrO2), yttrium oxide (Y2O3) hafnium dioxide (HfO2), and aluminum oxide (Al2O3).

It is contraindicated to install the implant in the presence of acute inflammatory or infectious processes in live tissue, unsuitable bone volume and/or quality, systemic complications or diseases such as; bone metabolism disorders, blood clotting disorders, unsuitable healing capacity, incomplete jawbone growth, uncooperative and not motivated patient, abuse of drugs or alcohol, psychosis, prolonged functional disorders which resist any treatment with medications, xerostomia, weakened immunological system, diseases which require the use of steroids, endocrine diseases, insufficient oral hygiene, and pregnancy.

**HANDLING**

Perform drilling of the bone bed using drills in good cutting conditions and rotation between 500 and 800 rpm for bone type III / IV, and rotation between 800 and 1200 rpm for bone type I / II, under abundant irrigation. Select the sequence of drills according to the selected implant. The insertion depth of the drills must comply with the planning of the implant final position. For the installation of Zirconia Implants, perform it manually with the Torque Wrench Driver ou perform it with surgical motor using the maximum torque of 35 N.cm and drilling rotation of 30 rpm. Neodent Drills/Bone Taps must be used in accordance with their own instructions for use. For more information, consult the Neodent catalog.

Sequence for Zirconia Implants:

Ø Implant (mm)	Bone Type	Initial Drill	Drill 2.0	Drill 3.5	Drill 3.75	Countersink Drill for Zirconia 3.75	Bone Tap for Zirconia 3.75	Drill 4.3	Countersink Drill for Zirconia 4.3	Bone Tap for Zirconia 4.3
4.3	I / II	*	✓	✓				✓	✓	✓
	III	*	✓	✓				✓	✓	
	IV	*	✓	✓				✓		

Obs: In order to prepare the surgical alveolus after extraction, use sequences of the drill used in type I bone. For the mandible, use a bone tap.

\* Optional

✓ Recommended sequence

**NOTES:**

The Initial Drill and Drills 2.0, 3.5, 3.75 and 4.3 indicated in the drilling sequence can be from GM or Alvim lines.


Drill the bone bed according to the drill sequence indicated for the planned implant. The Bone Tap is recommended to finish the preparation of the surgical cavity in type I / II bone and should be used according to the following guidelines:

**Step 1:** To start the insertion with the Bone Tap, use the Contra-Angle handpiece. Fit the Bone Tap to the Hex Connection for Contra-Angle and configure the surgical motor with a rotation speed of a maximum of 30 rpm and a maximum torque of 35 N.cm. Start the engine and insert the Bone Tap into the surgical cavity, maintaining the drilling shaft until stability is obtained and / or it reaches the maximum torque of 35 N.cm.

**Step 2:** Then, proceed with the installation with the Torque Wrench, using the Hexagonal Connection for Torque Wrench. To do this, fit the Bone Tap to the Hexagonal Connection for Torque Wrench and perform insertion movements, clockwise, lightly pressing the Torque Wrench. If you encounter a lot of resistance, it is recommended to apply contratorque to continue with the insertion. Proceed with the insertion until the Bone Tap for Zirconia Implant reaches the mark corresponding to the chosen implant. To completely remove the Bone Tap from the surgical cavity, reverse the direction of the Torque Wrench counterclockwise and carefully remove it. If performed differently, removing the instrument can impair the thread formation.

### INSTALLATION OF THE DENTAL IMPLANT

**WARNING:** The images are merely illustrative and do not represent the actual dimensions and features of the product.

- The cardboard and blister packagings must be opened, manually, without the use of sterile gloves.
- Break the seal of the cardboard packaging and remove the blister.
- Open the blister pack.
- Deposit the sterile flask over the surgical field. **NOTE:** The clear tube and implant must be handled with a sterile surgical glove, in a surgical environment.
- Hold the bottle using the non-dominant hand and take the lid off.
- The internal support containing the implant and transfer piece must come out attached to the lid. To do so, remove the lid and the clear tube's internal support in the axial direction without making any lateral movements.
- Using the non-dominant hand, press the sides of the internal support promoting a "pincer effect" and immobilizing the implant.
- Keep the support pressed and remove the lid.
- For installation, capture the implant transfer piece with the Hexagonal Connection, keeping it stable and slightly rotating the internal support, searching for the perfect fit between connection and transfer piece.
- Take the transfer-implant assembly to the surgical cavity.
- In the surgical motor, use maximum torque of 35 N.cm and rotation of 30 rpm.
- Complete the installation of the implant to the bone level with the Torque Wrench Driver with the Hexagonal Connection for Torque Wrench, short or long.
- At the end of the installation, make sure that one of the six hexagonal faces of the transfer piece, corresponding to the implant indexer, is facing the vestibular face.
- After the installation, remove the transfer piece from the implant.

• Load application according to torque is described in the following table:

Load Application	Min. Torque (N.cm)	Max. Torque (N.cm)
Delayed Loading*	10	60
Immediate Loading	35	60

\* Associated with the use of Cover Screw or Healing Abutment.

### PROSTHETIC SEQUENCE

In the immediate postoperative period, according to the chosen loading protocol, install the cover screw, healing abutment or prosthetic abutment, respecting their corresponding indications and limitations.

### TRACEABILITY LABEL

This product is supplied along with three labels that allow its traceability and must be attached to the following documents: • Medical record; • Collection tax document; • Document to be delivered to the patient (ask your advisor). The identification and traceability are performed through numeric codes REF and LOT.

### PRESENTATION AND STERILIZATION

This product is for single use and supplied sterile by the ethylene oxide method, being unitarily packed in a packaging that offers triple protection: carton box, blister and clear tube.

### PRE-OPERATORY AND PROSTHETIC PLANNING

The model, diameter, length, position, and quantity of implants must be selected for each clinical case, considering the anatomy, the region to be rehabilitated, bone quality and quantity, as well as available space. Execute the wax-up diagnostic of the clinical case, ensuring the appropriate alignment between implant(s), prosthetic abutment(s), and prosthesis(es).

### PRECAUTIONS

- Due to limited mouth opening in the posterior region, an assessment is required prior to installation of long implants in the region of premolars and molars. The installation of long implants in these regions may require them to be installed at an angle.
- Inadequate surgical and/or prosthetic planning can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of abutments and/or prosthetic screws.
- Be aware in cases of patients that present signs of allergy or hypersensitivity to chemical elements of the material: stainless steel.
- Observe the conditions of the intraoral tissue, bone quality and quantity of the bed receiving the implant, through radiographic and/or tomography examinations. Not performing the pre-surgical assessment may compromise the success of the procedure.
- As for the systemic aspect, consider the general health of the patient. In particular, one must be careful in cases of patients who have allergies to drugs, local or systemic factors that may interfere with the healing process of the bone tissues or soft tissues, or the process of osseointegration. For example, the bone already exposed to radiation in the head and neck area, diabetes mellitus, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, previous treatment of jaws pathologies and oral mucosa abnormalities.
- Treatment with bisphosphonates results in potential risk of peri-implant osteonecrosis.
- Avoid approaching the proximities of the mandibular nerve channel during bone bed surgical preparation. Nerve damage may result in anesthesia, paresthesia, and dysthesia.

- The material to be used during the procedure must be sterile. Never use potentially contaminated items, since contamination may lead to infections.
- This product is for single use and cannot be re-sterilized. Neodent does not accept any responsibility for re-sterilized implants, regardless of who has carried out or how it was performed. Cleaning, disinfection, and sterilization of this product by oneself may compromise essential material and design characteristics.
- Reuse of this product may cause: adverse biological effects of residual products, microorganisms, and/or substances resulting from previous uses and/or reprocessing; changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put at risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.
- Do not use the product if its packaging is damaged.
- Do not use the product with its validity expired.
- The sterilization is guaranteed only if the sterile barrier (full blister) is not damaged.
- This product must be used immediately after opening of its packaging, at the moment of surgery. If it is not used, discard it.
- Follow the rules of asepsis when removing the implant from the sterile packaging.
- The maximum installation torque suggested is of 60 N.cm. Insertion torque greater than the one recommended can make the system inoperative. Applying an improper torque could lead to implant breakage or bone necrosis.
- Before each procedure, make sure the pieces are properly seated. Ensure the perfect fit between implant and transfer piece.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Check passivity and perform occlusal and interproximal adjustment after installation of the prosthesis, avoiding impairment of the implant/prosthesis assembly.
- Do not continue to insert the implant in case the transfer piece breaks. Insisting on the installation may cause damages to the system and/or the patient. Use the appropriate connection to remove the implant.
- Before each procedure, check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is any damage, deleted markings, compromised sharpening, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the dental surgeon's responsibility to use the Neodent products according to the instructions for use. • In case there is an accident related to the device, inform the manufacturer using the communication channels available on the website.

### ADVERSE EFFECTS

The installation of dental implants, as well as any other surgical procedure, may cause a slight discomfort and localized edema. More persistent symptoms can occur such as: chronic pain related to the dental implant, permanent paresthesia, dysesthesia, loss of maxillary/mandibular bone edge, systemic or localized infection, oroantral or oronasal fistula, adjacent teeth unfavorably affected, irreversible damage to adjacent teeth, fracturing of the implant, jaw, bone or prosthesis, esthetic problems, injury of the nerves, exfoliation, hyperplasia.

Any failure in osseointegration and loss of the prosthesis during treatment may be caused by: inadequate osteotomy, infections, diseases or systemic problems, low quality or volume of the remaining bone, absence or failure of irrigation, use of nonspecific instruments and/or instruments without cutting power, insufficient oral hygiene, occlusal trauma, lack of prosthetic passivity and lack of specific training.

### MAGNETIC RESONANCE IMAGING (MRI) – SAFETY INFORMATION

The Neodent Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neodent Implant System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

### POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need for professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. Evaluate the need for follow-up with periodic clinical appointments and radiographic examinations. The professional in charge is responsible for providing these guidelines.

### STORAGE CONDITIONS

This product must be stored, in its original packaging, in a clean and dry location and protected from direct sunlight.

### DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them. Before discarding them into the environment, it is recommended to observe the current legislation and adhere to it.

### DATE OF EXPIRATION

Written on the label.

### SYMBOLS

The table describes the symbols that may be printed on the product labeling. Please, refer to the physical labeling to see which symbols are applicable to the product.

Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
	Quantity
	Consult instructions for use
	Reprocessing forbidden
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community
	Importer in the European Community
	Swiss Authorised representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
	Notification required by FDA for United States market
	CE Mark

Symbols	Description
NON-STERILE 	Non-sterile

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Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries.  
Please, contact the authorized distributor.

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© 2021 – JJGC Indústria e Comércio de Materiais Dentários S.A. All rights reserved.
© 2021 – JJGC Indústria e Comércio de Materiais Dentários S.A. Tous droits réservés.